

Stuart, SRT1 and SRT2 Roller Mixers

- // Rocking and rolling action for complete mixing
- // Ideal for mixing blood samples
- // Removable rollers to accommodate large bottles
- // Designed for continuous operation between 4°C and 60°C
- // Choice of two sizes, 5 rollers or 10 rollers
- // 3 year warranty



Stuart, SRT1 and SRT2 Roller Mixers

The roller mixers provide a gentle, but highly effective, rocking and rolling action. The rollers rotate at 33rpm and rise and fall 16mm.

Recommended for mixing blood samples, viscous substances and liquid-solid suspensions where minimum aeration is required. Both roller mixers are robustly constructed and designed for easy cleaning, having plastic rollers and a stainless steel drip tray to collect accidental spillages. Most sizes of tubes, Bijoux vials, Universals and bottles can be accommodated.

For large capacity bottles (up to 120mm diameter), one or more of the rollers is easily removed to give improved grip. Similarly, it is possible to quickly convert the units to provide a 'roll only' action and function as a simple bottle roller.

Specifications

Model	SRT1	SRT2
No. of rollers	5	10
Roller diameter/length, mm	30/340	30/340
Speed per minute, rpm	33	33
Orbit dia./amplitude, mm	16	16
Dimensions (wxdxh), mm	200x500x90	380x500x90
Weight, kg	5	7

Cat. No.	Description
364-061	SRT1 Roller mixer (5 roller model)
364-069	SRT2 Roller mixers (10 roller model)

364-069

