

Containment trays and APET anti-static liners

Products Company Contact Home Buy online



Containment Trays and APET Anti-Static Liners

Trays are rigid PVC with rounded corners for easy cleaning, they provide a safety containment area for working with hazardous substances. They are available in the radiation safety colour of yellow and display the radiation Trefoil or in white for use outside of radiation areas.

The liners are constructed of APET, which is environmentally friendly when incinerated. They are non-porous, anti-static and can be rolled up for disposal.

Radiation hazard containment trays

Cat. No.	Description	External dimensions, mm	Internal dimensions, mm	Pack
526-139	Yellow tray, base only	460x260	449x249	1
526-140	Yellow tray, base only	540x340	529x329	1
526-141	Yellow tray, base only	570x540	559x529	1
526-142	Yellow tray, base only	680x540	669x529	1
526-143	Yellow tray, base only	700x460	689x449	1
536-144	Yellow tray, base only	1130x540	1119x529	1

General purpose containment trays

Cat. No.	Description	External dimensions, mm	Internal dimensions, mm	Pack
526-145	White tray, base only	460x260	449x249	1
526-146	White tray, base only	540x340	529x329	1
526-147	White tray, base only	570x540	559x529	1
526-148	White tray, base only	680x540	669x529	1
526-149	White tray, base only	700x460	689x449	1
526-150	White tray, base only	1130x540	1119x529	1

APET liners

Cat. No.	Description	External dimensions, mm	Internal dimensions, mm	Pack
526-133	Anti-static liners	460x260	449x249	25
526-134	Anti-static liners	540x340	529x329	25
526-135	Anti-static liners	570x540	559x529	25
526-136	Anti-static liners	680x540	669x529	25
526-137	Anti-static liners	700x460	689x449	25
526-138	Anti-static liners	1130x540	1119x529	25

Base trays suitable for biohazard work are also available, please contact our Technical Office on 01342 826836 for details.