

B-Safe Bin Box Bags in Clear Un-badged Polyethylene

Products

Company

Contact

Home

Buy online

526-089 to 526-092



B-Safe Bin Box Bags in Clear Un-badged Polyethylene

Supplied in benchtop dispenser boxes. Heat sealed for leak protection. Each bag is box shaped and has a double string neck pull. No printing. Polyethylene (500 gauge).

Cat. No.	Size (WxDxH)	Pack size
526-089	150x150x150mm	25
526-090	200x200x250mm	25
526-091	280x280x510mm	25
526-092	490x270x400mm	25