

Cooler Ice Buckets

[Products](#)

[Company](#)

[Contact](#)

[Home](#)

[Buy online](#)

190-101



Cooler Ice Bucket

The Cooler is specially designed to provide a rugged container for efficient thermal storage of ice and dry ice on the laboratory bench. The smooth, easy-to-clean, chemically resistant, white, high impact, polystyrene body and lid have a 25mm thick double insulated wall filled with expanded polystyrene foam. Each Cooler is made with four convenient support holes to accept your single channel pipettors. Ideal for laboratory use, the Cooler Ice Bucket is stable, leakproof and has a capacity of 3.5 litres.

External dimensions

280 x 300 x 200 (hwxwd), mm

Cat No.	Description
190-101	Cooler ice bucket
190-102	Lid for Cooler ice bucket