

Safety Zippette™, bottle-top dispenser

Products

Company

Contact

Home

Buy online



268-005

The MicroZippette™, Hand-held Dispenser

The MicroZippette™ is a spring-loaded, self-priming, hand-held dispenser. Ideal for rapid multi-dispensing of reagents as an alternative to micropipettors. Four dual fixed volume models are available: 100/250µl, 500/750µl, 500/1000µl and 1000/2000µl.

- Inaccuracy $\pm 1\%$, Imprecision 0.2%
- Positive lock volume adjustment
- Excellent precision and accuracy
- May be conveniently stored in a refrigerator
- Chemically resistant fluid pathway
- Supplied complete with dispense tube, intake tube and plastic coated safety bottle
- Fully autoclavable

Cat. No.	Description
268-005	MicroZippette™ 100 and 250µl volumes supplied with bottle
268-006	MicroZippette™ 500 and 750µl volumes supplied with bottle
268-029	MicroZippette™ 500 and 1000µl volumes supplied with bottle
268-007	MicroZippette™ 1000 and 2000µl volumes supplied with bottle