

Safetypette™, pipette filler

Products

Company

Contact

Home

Buy online



Safetypette, Pipette Filler

A simple yet controllable pipette filler, the Safetypette is designed to give superior pipetting control to left and right-handed users. The lightweight Safetypette sits comfortably in the hand, enabling smooth control in both aspirating and dispensing, using 1ml to 100ml pipettes. Simply squeeze the silicone bulb, and the thumb lever controls both the aspirate and dispense speed. Incorporates central button to blow out residual liquid in blow-out pipette. Supplied with three additional coloured pipette cones.

- Comfortable and simple to use
- Precise pipetting control
- Robust and lightweight
- Large capacity bulb
- Compatible with blow-out pipettes
- Simple thumb button operation
- Uses replaceable standard membrane filter
- Autoclavable silicone pipette holder
- Supplied with integral 0.45µm filter and coloured pipette cones for safety colour coding in white, mint, mauve and pink
- Easily maintained and cleaned

Cat. No.	Description	Pack
475-102	Safetypette pipette filler with integral 0.45µm filter and coloured pipette cones	–
266-004	Replacement silicone pipette holder	1
266-003	Replacement filter set 0.45µm	5
266-089	Replacement filter set 0.2µm	5
266-006	Replacement sterile filter 0.45µm	1